

CGA V-5—2008

REAFFIRMED 2013

**DIAMETER INDEX SAFETY SYSTEM
(NONINTERCHANGEABLE LOW
PRESSURE CONNECTIONS FOR
MEDICAL GAS APPLICATIONS)**

SIXTH EDITION



**COMPRESSED GAS
ASSOCIATION, INC.**

PLEASE NOTE:

The information contained in this document was obtained from sources believed to be reliable and is based on technical information and experience currently available from members of the Compressed Gas Association, Inc. and others. However, the Association or its members, jointly or severally, make no guarantee of the results and assume no liability or responsibility in connection with the information or suggestions herein contained. Moreover, it should not be assumed that every acceptable commodity grade, test or safety procedure or method, precaution, equipment or device is contained within, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedure.

This document is subject to periodic review, and users are cautioned to obtain the latest edition. The Association invites comments and suggestions for consideration. In connection with such review, any such comments or suggestions will be fully reviewed by the Association after giving the party, upon request, a reasonable opportunity to be heard. Proposed changes may be submitted via the Internet at our web site, www.cganet.com.

This document should not be confused with federal, state, provincial, or municipal specifications or regulations; insurance requirements; or national safety codes. While the Association recommends reference to or use of this document by government agencies and others, this document is purely voluntary and not binding unless adopted by reference in regulations.

A listing of all publications, audiovisual programs, safety and technical bulletins, and safety posters is available via the Internet at our website at www.cganet.com. For more information contact CGA at Phone: 703-788-2700, ext. 799. E-mail: customerservice@cganet.com.

Work Item 12-055
Medical Equipment Committee

NOTE—Technical changes from the previous edition are underlined.

NOTE—Appendices A, B, C, D, E, F, G and H (Normative) are a requirement.

NOTE—No technical information has been changed from the 2008 edition. This reaffirmed edition may include minor editorial changes.

REAFFIRMED: 2013
SIXTH EDITION: 2008
FIFTH EDITION: 2005
FOURTH EDITION: 2000

© 2008 The Compressed Gas Association, Inc. All rights reserved.

All materials contained in this work are protected by United States and international copyright laws. No part of this work may be reproduced or transmitted in any form or by any means, electronic or mechanical including photocopying, recording, or any information storage and retrieval system without permission in writing from The Compressed Gas Association, Inc. All requests for permission to reproduce material from this work should be directed to The Compressed Gas Association, Inc., 14501 George Carter Way, Suite 103, Chantilly VA 20151. You may not alter or remove any trademark, copyright, or other notice from this work.

Contents	Page
1 Introduction.....	1
2 Scope	1
3 CGA V-5 revisions	1
3.1 Second edition (1978)	1
3.2 Third edition (1989)	1
3.3 Fourth edition (2000).....	2
3.4 Fifth edition (2005).....	2
3.5 Sixth edition (2008).....	2
4 Description of the Diameter Index Safety System	2
4.1 General features.....	2
4.2 Optional features	3
5 Connection assignments	3
5.1 Connection Nos. 1020 through 1200 (1000 series).....	3
5.2 Connection Nos. CGA 1020-A through 1200-A (1000-A series)	4
5.3 Connection No. CGA 1220 for suction (vacuum)	4
5.4 Connection No. CGA 1240 for oxygen	5
5.5 Connection Nos. CGA 2020 through 2160 (2000 series).....	5
5.6 Connection No. CGA 2220 for waste anesthetic gas disposal (gas evacuation))	5
6 References	7

Table

Table 1—Medical gas connection numbers assigned to gases and gas mixtures in the Diameter Index Safety System	6
---	---

Figures

Figure 1—Typical standard connection, Series No. 1000.....	4
Figure 2—Typical co-standard connection, Series No. 1000-A	4

Appendices

Appendix A—Drawing No. 1000 (Normative).....	8
Appendix B—Drawing No. 1000-A (Normative)	9
Appendix C—Connection Nos. 1020 through 1200-A (Normative)	10
Appendix D—Connection No. CGA 1220 (Normative).....	11
Appendix E—Connection No. CGA 1240 (Normative).....	12
Appendix F—Drawing No. 2000 (Normative).....	13
Appendix G—Connection Nos. CGA 2020 through 2160 (Normative).....	14
Appendix H—Connection No. CGA 2220 (Normative).....	15

This page is intentionally blank.

1 Introduction

The administration of medical gases to patients is subject to a variety of hazards that can result from human error due to preoccupation, a mental lapse, or carelessness. This publication was developed to reduce the possibility of inadvertent substitution of the wrong medical gases by personnel when using equipment for administering anesthesia, resuscitation, and therapy.

CGA V-1, *Standard for Compressed Gas Cylinder Valve Outlet and Inlet Connections*, establishes assigned standard connections for threaded outlets on cylinder valves, which results in the removal from service of all cylinder valves having nonconforming threaded outlets [1].¹ CGA V-1 also describes the details of the Pin Index Safety System for flush-type medical gas cylinder valves. This system reduces the possibility of inadvertent substitution of the wrong gases on equipment with yoke-type connections. The Pin Index Safety System has been adopted by national and international standards organizations including the Canadian Standards Association International (CSA) and the International Organization for Standardization (ISO).

2 Scope

The Diameter Index Safety System defines noninterchangeable removable connections for use with medical gases at pressures of 200 psig (1380 kPa) or less.² Removable, threaded connections are readily and commonly engaged or disengaged in routine use and service and are often found on regulators, resuscitation equipment, anesthesia equipment, therapy equipment, and other similar equipment.

The Diameter Index Safety System supplements but does *not* replace:

- any of the means for medical gas identification currently in use;
- the Pin Index Safety System;
- the existing threaded outlet standards for cylinder valves; or
- automatic quick coupler valves that also provide noninterchangeable connections for medical gases and suction equipment.

CGA does not designate where the Diameter Index Safety System should be applied.

3 CGA V-5 revisions

3.1 Second edition (1978)

The second edition of CGA V-5 included the original twelve connections described in the 1959 edition and added a co-standard 1000-A series. Eight new connections (1500 series) designed and shown in SI dimensions were also added.

3.2 Third edition (1989)

3.2.1 1500 series replaced by 2000 series

The third edition abandoned the 1500 series connections in favor of a new smaller 2000 series, which offered eight new connections that are designed and shown in customary U.S. dimensions (inches).

3.2.2 Revised nitrous oxide-oxygen mixture connection

Connection No. CGA 1570, which was described in the 1978 edition and assigned to nitrous oxide-oxygen mixtures (nitrous oxide 47.5% to 52.5%), was replaced by Connection No. CGA 2020.

¹ References are shown by bracketed numbers and are listed in order of appearance in the reference section.

² kPa shall indicate gauge pressure unless otherwise noted as (kPa, abs) for absolute pressure or (kPa, differential) for differential pressure. All kPa values are rounded off per CGA P-11, *Metric Practice Guide for the Compressed Gas Industry* [2].