



**CGA V-7.1—2021  
STANDARD METHOD OF  
DETERMINING CYLINDER  
VALVE OUTLET  
CONNECTIONS FOR  
MEDICAL GASES**

**THIRD EDITION**

**PLEASE NOTE:**

The information contained in this document was obtained from sources believed to be reliable and is based on technical information and experience currently available from members of the Compressed Gas Association, Inc. and others. However, the Association or its members, jointly or severally, make no guarantee of the results and assume no liability or responsibility in connection with the information or suggestions herein contained. Moreover, it should not be assumed that every acceptable commodity grade, test or safety procedure or method, precaution, equipment or device is contained within, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedure.

This document is subject to periodic review, and users are cautioned to obtain the latest edition. The Association invites comments and suggestions for consideration. In connection with such review, any such comments or suggestions will be fully reviewed by the Association after giving the party, upon request, a reasonable opportunity to be heard. Proposed changes may be submitted via the Internet at our website, [www.cganet.com](http://www.cganet.com).

This document should not be confused with federal, state, provincial, or municipal specifications or regulations; insurance requirements; or national safety codes. While the Association recommends reference to or use of this document by government agencies and others, this document is purely voluntary and not binding unless adopted by reference in regulations.

A listing of all publications, audiovisual programs, safety and technical bulletins, and safety posters is available via the Internet at our website at [www.cganet.com](http://www.cganet.com). For more information contact CGA at Phone: 703-788-2700, ext. 799. E-mail: [customerservice@cganet.com](mailto:customerservice@cganet.com).

Work Item 16-014  
Medical Gases Committee

---

NOTE—Technical changes from the previous edition are underlined.

NOTE—Appendices A, B, and C (Informative) are for information only.

THIRD EDITION: 2021  
SECOND EDITION: 2011  
REAFFIRMED: 2008  
REAFFIRMED: 2003  
FIRST EDITION: 1997

© 2021 The Compressed Gas Association, Inc. All rights reserved.

All materials contained in this work are protected by United States and international copyright laws. No part of this work may be reproduced or transmitted in any form or by any means, electronic or mechanical including photocopying, recording, or any information storage and retrieval system without permission in writing from The Compressed Gas Association, Inc. All requests for permission to reproduce material from this work should be directed to The Compressed Gas Association, Inc., 8484 Westpark Drive, Suite 220, McLean VA 22102. You may not alter or remove any trademark, copyright or other notice from this work.

<b>Contents</b>	<b>Page</b>
1 Introduction.....	1
2 Scope and Purpose.....	1
2.1 Scope .....	1
2.2 Purpose .....	1
3 Definitions.....	1
4 Description of the method .....	2
5 Principles for connection assignment selection.....	2
5.1 Single medical gases assigned connections in CGA V-1.....	2
5.2 Single medical gases not assigned a connection in CGA V-1 or medical gas mixtures.....	3
6 References .....	6
<b>Appendices</b>	
Appendix A—Selection flow diagram (Informative).....	7
Appendix B—Medical gas/gas mixture algorithm for container valve outlet and inlet connection assignment (Informative).....	10
Appendix C—Connections assigned by CGA V-1, ISO 407 and derived from CGA V-7.1 (Informative).....	14

This page is intentionally blank.

## 1 Introduction

CGA V-1, *Standard for Compressed Gas Cylinder Valve Outlet and Inlet Connections* includes valve outlet connections for more than 175 different pure gases and chemicals but includes only a few gas mixtures that have a history of commercial or medical importance for which standards were in demand [1].<sup>1</sup>

In response to the increase in the number of industrial gas mixtures, CGA published CGA V-7, *Standard Method of Determining Cylinder Valve Outlet Connections for Industrial Gas Mixtures* [2]. CGA V-7 excludes medical gas mixtures.

Medical gases are given special treatment because of their applications involving medical patients. CGA determined that a method of selecting cylinder valve outlet connections for medical gases would help manufacturers of medical gases and medical equipment select appropriate connections for new medical gases and use existing connection assignments better. This method addresses known potential hazards to the user or the patient caused by inappropriate connections. If a medical gas equipment or medical gas manufacturer believes the existing connection assignment does not adequately address safety requirements, that manufacturer shall request a new connection assignment from CGA as detailed in CGA V-1 [1].

Where used in this standard, the term medical gases shall mean compressed gases in cylinders that are labeled as either drugs or medical devices, including gas mixtures. The definition of drugs and medical devices as it relates to these gases shall be that as promulgated by the Food and Drug Administration (FDA) in the United States and by the Health Canada Therapeutic Products Programme in Canada.

## 2 Scope and Purpose

### 2.1 Scope

This standard applies to the selection of an appropriate cylinder valve outlet connection for medical gases for pressures up to 3000 psi (20 680 kPa) at 70 °F (21.1 °C).<sup>2</sup> It does not apply to industrial gas mixtures, cylinders manifolded in transport and use, or to mixtures at higher pressures. Both threaded and pin-indexed connection assignments are included with reference to V-1 [1].

This revision replaces the use of limited standards with alternate standards for the purpose of assuring only one connection assignment will be appropriate. This revision brings consistency with ISO 407:2004, *Small medical gas cylinders—Pin-index yoke-type valve connections*, to standardize pin-indexed yoke-type valve connections globally [4].

### 2.2 Purpose

The purpose of this standard is to establish a procedure for selecting the appropriate cylinder valve outlet connection for medical gases. It does not attempt to identify which medical gases may be safely manufactured nor indicate that the connection is safe or effective for any medical procedure. It is the responsibility of the manufacturer, the medical community, and appropriate regulatory bodies to make those decisions.

## 3 Definitions

For the purpose of this publication, the following definitions apply.

### 3.1 Publication terminology

#### 3.1.1 Shall

Indicates that the procedure is mandatory. It is used wherever the criterion for conformance to specific recommendations allows no deviation.

---

<sup>1</sup> References are shown by bracketed numbers and are listed in order of appearance in the reference section.

<sup>2</sup> kPa shall indicate gauge pressure unless otherwise noted as (kPa, abs) for absolute pressure or (kPa, differential) for differential pressure. All kPa values are rounded off per CGA P-11, *Guideline for Metric Practice in the Compressed Gas Industry* [3].